

## Trans-Femoral Prosthesis (Above Knee Prosthesis)



Your prosthesis has be created specifically for you. Always follow specific instructions given to you by your doctor and practitioner.

## Wearing & Break-in Schedule

Always check your skin after wearing your device. Any areas of pink should fade away after 10-15 minutes. If marks do not disappear after 10-15 minutes, or if you observe deep red marks, blisters, bruises, or sores, contact your practitioner immediately. If you have any concerns **do not** wear the device until you follow up with your practitioner.

It is important for you to follow the "break-in" schedule of wearing. Following this schedule allows time for your tissues to adjust to the new pressures on your bone and skin.

<u>Break-in Schedule:</u> Begin wearing your prosthesis for 15-30 minute intervals. Remove it earlier if you notice any pain or discomfort. Remove the prosthesis and check your skin for any redness that does not go away within 10-15 minutes. If marks do not disappear after 10-15 minutes, or if you observe deep red marks, blisters, bruises, or sores, contact your practitioner immediately. Gradually increase the time worn until comfortable full time wear is achieved using the following schedule:

Day 1 – Wear for 1 hr.

Day 2 – Wear for 2 hrs.

Day 3 – Wear for 3 hrs.

Day 4 – Wear for 4 hrs.

Day 5 – Wear for 5 hrs.

Day 6 – Wear for 6 hrs.

Day 7 – Wear for 7 hrs.

Day 8 – Wear for 8 hrs.

Day 9 – Wear for 9 hrs.

## How to put on your prosthesis

Prior to wearing your prosthesis, it is important to check the following things on your device and residual limb.

- Make sure your skin is <u>clean and dry</u>.
- Avoid the use of skin creams and salves prior to wearing the device.
- Examine your skin and look for abnormal redness, blisters, bruises, or sores. Using a mirror may be helpful if you have difficulty seeing any areas on the back of your leg.
- Make sure there are no cracks, rough edges, or worn areas on the device.
- Make sure any straps are secure on the device.
- Check your prosthetic socks and liner (if applicable), for worn areas or holes.

<sup>\*\*</sup> Continue to add one hour each day until full wear is achieved. This usually takes at least two weeks.